

EC View

ECG editing station



CARDIOLEX

EASY TO USE ECG

EC View

(Catalogue no: CN V-A)

Streamline
your digital
workflow!

EC View is a modern editing station and decision management system for the diagnosis of Rest ECG, Exercise ECG and Long ECG. It is used for viewing/editing examinations that are stored on a central storage system. EC View can be used by multiple license users simultaneously and the same ECG examination can be viewed at several PCs.

BENEFITS WITH EC VIEW



OVERVIEW



EASY TO USE



TIME MANAGEMENT



TEMPLATES




FULL DISCLOSURE

OVERVIEW


Using combined search and sorting, the operator can easily find the examination results that are relevant. Once an examination is opened the ECG is presented in different tabs, providing an overview of the ECG. Buttons for save, print and back to search are visible in all tabs. The different tabs are:

TABS	Rest ECG	Exercise ECG	Long ECG
Diagnosis	✓		
Beats	✓		
Leads	✓		
Serial report	✓		✓
Trend	✓	✓	✓
Report	✓	✓	
Summary		✓	✓
Rhythm		✓	
Arythmia/Event		✓	
Complex		✓	
ECG			✓



EASY TO USE

EC View uses authorization based on active directory with single-sign-on, resulting in fast and secure access for the operator. The recognition factor for the operator is high since the presentation of ECG is identical regardless of the source, streamlining the assessments. EC View holds bidirectional communication for retrieval of ECG examinations from EC Store. The simplified workflow decreases examination time and prevents manual entry mistakes.




TIME MANAGEMENT

The longer examinations, Exercise ECG and Long ECG, have a time reference marker. It is a useful tool for searching and reviewing events on the different pages and tabs. The time reference is synchronized for the same time on all trends: by setting the time reference in one location in the ECG it is possible to review other parameters, complexes and leads for the same interval of time. Events can be marked and are presented in the report.



TEMPLATES

By using report templates defined by the user, the workflow efficiency for the operator increases. Report templates can be predefined in Settings or set temporarily by the operator using a pop-up menu. It is for example possible to preconfigure test as templates to make it easier to enter commonly used text in comment boxes.



FULL DISCLOSURE

Complete records of the examinations can be viewed in high quality so that no information is lost. Patient information, reports, critical values, notes etc. are stored along with the examination, providing full disclosure.



CARDIOLEX

EASY TO USE ECG

Cardiolex Medical AB develops and markets easy to use systems for ECG.

Our ECG management system and various ECG systems exist on 2/3 of Swedish hospitals, which means that the system has over 10 000 users in Sweden. We offer the complete solution within ECG, from the systems receiving the ECG signal from the patient to the ECG management systems, storing and managing generated data. Our products work for a wide range of users, from University hospitals to small clinics. Our focus is on creating ECG solutions that are user-friendly, reliable and secure.

Read more about our products on www.cardiolex.com.



Cardiolex Medical AB
Vretenvägen 6, SE-171 54 Solna
+46 (0) 8 564 820 50 · info@cardiolex.se
www.cardiolex.com